

Product data sheet

GENERAL INFORMATION

Product: Centrifuge tube, polypropylene, with graduation, 50 ml

Description: Conical bottom without skirt. Autoclavable. Dimensions: tube \varnothing 29 mm, 115 mm high and cap \varnothing 33,5 mm. Withstands RCF of 10000xg. Successfully passed 95 kPa pressure test.

SPECIFICATIONS

Reference	pcs/ box	capacity	type
CTGP-050-050	50	50 ml	non-sterile
"Saving package" CTGP-050-500	500	50 ml	non-sterile

PACKING

Type: Carboard box

Label:

labbox	CTGP-050-050
	Graduated Centrifuge tubes with screw cap, 50 ml, PP, 50 pcs
	Tubos de centrifuga graduados con tapón de rosca, 50 ml, PP, 50 uds
	Tubes à centrifuger gradués en polypropylène, 50 ml, non stériles, 50 unit.
	Batch n ^o .: xxxx

MATERIAL

Polypropylene

Good chemical resistance and autoclavable at 121°C.

