

S.C.A.T. Europe Sets Overview

HPLC Supply and Waste

Generation V2.0



när special
är standard



THE System for
Solvent Safety

HPLC Starter Kits: Supply

- with 4x SafetyCaps incl.:
- 4x Air valves
 - Fittings
 - Blind plugs



HPLC Supply-Sets: Supply

- comprising:
- SafetyCap incl.: Air valve, Fittings, Blind plugs
 - Laboratory glass bottle
 - Capillary / Capillaries
 - Suction filters



HPLC Waste-Sets: Waste

- comprising:
- SafetyWasteCap incl.: Fittings, Tube connector, Blind plugs
 - Exhaust filter
 - Canister



HPLC Supply and Waste Set: Supply and Waste - Complete Set - 1 Set - suitable for all HPLC systems!

- comprising:
- 4x SafetyCaps incl. Air valves, Fittings, Blind plugs
 - 4x Laboratory glass bottles
 - 1x SafetyWasteCap incl. Fittings, Tube connectors, Blind plugs
 - 1x Exhaust filter
 - 1x Canister



Part No.	Description
399200	HPLC Starter Kit 1 - Supply , V2.0: 3x SafetyCaps I (307019), 1x SafetyCap II (307909), 1x blind plug, 5x fittings 1.6/2.3/3.2 mm, 4x air valves (317010)
399201	HPLC Starter Kit 2 - Supply , V2.0: 4x SafetyCaps II (307909), 4x blind plugs, 8x fittings 1.6/2.3/3.2 mm, 4x air valves (317010)
399202	HPLC Starter Kit 3 - Supply , V2.0: 4x SafetyCaps III (307910), 8x blind plugs, 12x fittings 1.6/2.3/3.2 mm, 4x air valves (317010)
307300	HPLC Supply-Set I , V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle round & clear glass + 1x 1.5 m capillary 3.2 mm OD, 1x PTFE suction filter + air valve (317010)
307301	HPLC Supply-Set I , V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle square & clear glass + 1x 1.5 m capillary 3.2 mm OD, 1x PTFE suction filter + air valve (317010)
307312	HPLC Supply-Set I , V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle round & brown glass + 1x 1.5 m capillary 3.2 mm OD, 1x PTFE suction filter + air valve (317010)
307303	HPLC Supply-Set II , V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle round & clear glass + 2x 1.5 m capillaries 3.2 mm OD, 2x PTFE suction filters + air valve (317010)
307305	HPLC Supply-Set II , V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle square & clear glass + 2x 1.5 m capillaries 3.2 mm OD, 2x PTFE suction filters + air valve (317010)
307313	HPLC Supply-Set II , V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle round & brown glass + 2x 1.5 m capillaries 3.2 mm OD, 2x PTFE suction filters + air valve (317010)
307304	HPLC Supply-Set III , V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle round & clear glass + 3x 1.5 m capillaries 3.2 mm OD, 3x PTFE suction filters + air valve (317010)
307306	HPLC Supply-Set III , V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle square & clear glass + 3x 1.5 m capillaries 3.2 mm OD, 3x PTFE suction filters + air valve (317010)
307314	HPLC Supply-Set III , V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle round & brown glass + 3x 1.5 m capillaries 3.2 mm OD, 3x PTFE suction filters + air valve (317010)
307347	HPLC Supply-Set III , V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle Youtility & clear glass + 3x 1.5 m capillaries 3.2 mm OD, 3x PTFE suction filters + air valve (317010)
307307	HPLC Waste-Set 1 , V2.0: SafetyWasteCap GL45, 2x 2.3/3.2 mm OD, 1x 6.4 - 9 mm ID (307923), 5 liter canister, PE-HD (107951), exhaust filter M, with splash guard and exchange indicator (310535)
307310	HPLC Waste-Set 2 , V2.0: SafetyWasteCap GL45, 2x 2.3/3.2 mm OD, 1x 6.4 - 9 mm ID (307923), 10 liter canister, PE-HD (107952), exhaust filter M, with splash guard and exchange indicator (310535)
307328	HPLC Waste-Set 3 , V2.0: SafetyWasteCap S50, 2x 3.2 mm OD, 1x 5.0 - 11.5mm ID (108025), 5 liter space saving canister, PP (107998), exhaust filter S, with splash guard and exchange indicator (310534)
307337	HPLC Supply and Waste Set, V2.0, suitable for all HPLC systems